λ polyA⁺ RNA-A

Code No.TX802

Size : 10 μ l

Shipping at — 80°C

Store at — 80°C

Supplied Reagents : DEPC treated water 1 ml

Lot No.

Conc. : 10 ng/ μ l

Expiration Date :

Description:

λ polyA⁺ RNA-A is RNA fragments (about 1.0 kb) attached with poly A, which are synthesized by *in vitro* transcription using λ DNAs template. These RNA fragments have the sequence complementary with λ DNA fragments (λ 1000 or lambda A) which are spotted on the products of TaKaRa IntelliGene™ Series. Store this product separately in each small amount after diluting them to the appropriate concentration by DEPC treated water.

For normalization of DNA microarray for eukaryote, this product is utilized for internal control.

For Real-time RT-PCR for eukaryote, this product is utilized for external RNA as reference with the use of Real Time Primer for λ poly A (Cat. #3786)

Usage :

Hybridization of DNA microarray

Before the fluorescent-labeling reaction, add the appropriate amount* of λ polyA⁺ RNA-A into RNA sample. Then label the sample and use it as the probe for hybridization. Adding the same amount of λ poly A⁺ RNA-A into both reference and test sample will enable to normalize the signals after the labeling reaction. Store this product in aliquots after diluting them to the appropriate concentration by DEPC treated water.

* : In case of using the TaKaRa IntelliGene™ for Eukaryote, it is recommended to add 100 - 500 pg of λ poly A⁺ RNA-A into each RNA sample.

To utilize this product for self-production of DNA microarray, immobilize DNA fragments (λ 1000 or λ DNA-A), which are produced by using λ Control Template & Primer Set-A (Cat. #TX803), on the DNA chips.

Real time RT-PCR

Regarding protocol, please refer to the product manual supplied with Real Time Primer for λ poly A (Cat.# 3786). Dilute this product by Easy Dilution (supplied in Cat.#3786) as protocol.

Notice :

Avoid freeze-thaw of this product , which are liable to disintegrate. Utilize this product according to the Notice for using RNA samples.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com.

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